

Datasheet

PRMT5 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00010419-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human PRMT5 protein.

Immunogen: PRMT5 (NP_006100.2, 1 a.a. ~ 637 a.a) full-length human protein.

Sequence:

MAAMAVGGAGGSRVSSGRDLNCVPEIADTLGAVAKQ
GFDFLCMPVFHPRFKREFIQEPAKNRPGPQTRSDLLL
SGRDWNTLIVGKLSPWIRPDSKVEKIRRNSEAAMLQEL
NFGAYLGLPAFLPLNQEDNTNLARVLTNHIHTGHHSS
MFWMRVPLVAPEDLRDDIIENAPTTHTTEYSGEEKTW
MWWHNFRTLCDYSKRIAAVALEIGADLPSNHVIDRWLG
EPIKAAILPTSIFLTNKKGFPVLSKMHQRLIFRLLKLEVQ
FIITGTNNHSEKEFCSYLQYLEYLSQNRPPPNAYELFA
KGYEDYLQSPLQPLMDNLESQTYEVFEKDPIKYSQYQ
QAIYKCLLDREVPEEEKDTNVQVLMVLGAGRGPLVNAS
LRAAKQADRRIKLYAVEKNPNAVVTLENWQFEEWGS
QVTVVSSDMREWVAPEKADIIVSELLGSFADNELSPEC
LDGAQHFLKDDGVSIPEYTSFLAPISSSKLYNEVRAC
REKDRDPEAQFEMPYVVRHLNHFHQLSAPQPCFTFSH
PNRDPMDNNRYCTLEFPVEVNTVLHGFAGYFETVLY
QDITLSIRPETHSPGMFSWFPILFPIKQPITVREGQTICV
RFWRCSNSKKVWYEWAVTAPVCSAIHNPTGRSYTIGL

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 10419

Gene Symbol: PRMT5

Gene Alias: HRMT1L5, IBP72, JBP1, SKB1, SKB1Hs