

Datasheet

MFGE8 purified MaxPab rabbit polyclonal antibody (D01P)

Catalog Number: H00004240-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human MFGE8 protein.

Immunogen: MFGE8 (NP_005919.1, 1 a.a. ~ 387 a.a) full-length human protein.

Sequence:

MPRPRLAALCGALLCAPSLLVALDICKSNPCHNGGLC
EEISQEVRGDVFPSTCTCLKGYAGNHCETKCVEPLG
MENGNIANSQIAASSVRVTFGLQHWVPELARLNLAG
MVNAWTPSSNDDNPWQVNLRRMWVTGVVTQGAS
RLASHEYLKAFKVAYSLNGHEFDHVDNKKHKEFVGN
WNKNVHVNLFETPVEAQYVRLYPTSCHTACTLRFEL
LGCELNGCANPLGLKNNSIPDKQITASSSYKTWGLHLF
SWNPSYARLDKQGNFNAWVAGSYGNDQWLQVDLGS
SKEVTGIITQGARNFGSVQFVASYKVAYSNDSANWTE
YQDPRTGSSKIFPGNWDNHSKKLNFETPILARYVRIL
PVAWHNRALRLELLGC

Host: Rabbit

Reactivity: Human

Applications: WB-Ti, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 4240

Gene Symbol: MFGE8

Gene Alias: BA46, EDIL1, HsT19888, OAcGD3S, SED1, SPAG10, hP47