

Datasheet

MGMT MaxPab rabbit polyclonal antibody (D01)

Catalog Number: H00004255-D01

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human MGMT protein.

Immunogen: MGMT (NP_002403.1, 1 a.a. ~ 207 a.a) full-length human protein.

Sequence:

MDKDCCEMKRTTLDSP LGKLELSGCEQGLHEIKLLGKG
TSAADAVEVPAPAAVLGGPEPLMQCTAWLNAYFHQP
EAIEEFPVPALHHPVFQQESFTRQVLWKLLKVVKFGEV
ISYQQLAALAGNPKAARAVGGAMRGNPVPILIPCHRVV
CSSGAVGNYSGLAVKEWLLAHEGHRLGKPG LGGSS
GLAGAWLKGAGATSGSP PAGRN

Host: Rabbit

Reactivity: Human, Mouse

Applications: IP, WB-Ti, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 4255

Gene Symbol: MGMT

Gene Alias: -