

## Datasheet

### MT2A monoclonal antibody (M01), clone 6G2

**Catalog Number:** H00004502-M01

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody  
raised against a full-length recombinant MT2A.

**Clone Name:** 6G2

**Immunogen:** MT2A (NP\_005944.1, 1 a.a. ~ 61 a.a)  
full-length recombinant protein with GST tag. MW of the  
GST tag alone is 26 KDa.

**Sequence:**

MDPNCSCAAGDSCTCAGSCKCKECKCTSCKKSCCSC  
CPVGCAKCAQGCICKGASDKCSCCA

**Host:** Mouse

**Reactivity:** Human

**Applications:** ELISA, IF, S-ELISA, WB-Re  
(See our web site product page for detailed applications  
information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product  
page for detailed protocols

**Isotype:** IgG2a Lambda

**Storage Buffer:** In 1x PBS, pH 7.4

**Storage Instruction:** Store at -20°C or lower. Aliquot to  
avoid repeated freezing and thawing.

**Entrez GeneID:** 4502

**Gene Symbol:** MT2A

**Gene Alias:** MT2