

Datasheet

NAP1L3 purified MaxPab rabbit polyclonal antibody (D01P)

Catalog Number: H00004675-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human NAP1L3 protein.

Immunogen: NAP1L3 (AAH34954.1, 1 a.a. ~ 506 a.a) full-length human protein.

Sequence:

MAEADFKMVSEPVAHGVAEEEMASSTSDSGEESDSS
SSSSSTSDSSSSSTSGSSSGSGSSSSSGSTSSRSR
LYRKRVPEPSRRARRAPLGTNFVDRLPQAVRNRVQ
ALRNIQDECDKVDTLFLKAIHDLERKYAELNKPLYDRR
FQIINAEYEPTEECEWNSEDEEFSSDEEVQDNTPE
MPPLEGEEEEENPKENPEVKAEKEVPKEIPEVKDEEK
EVAKEITEVKAEEKADSKDCMEATPEVKEDPKEVPQV
KADDKEQPKATEAKARAAVRETHKRVPEERLQDSVDL
KRARKGKPKREDPKGIPDYWLIVLKNVDKLGPMIQKYD
EPILKFLSDVSLKFSKPGQPVSYTFEFHFLPNPYFRNE
VLVKTYIIKAKPDHNDPFFSWGWEIEDCKGCKIDWRR
GKDVTVTTTQSRTTATGEIEIQPRVVPNASFFNFFSPP
EIPMIGKLEPREDAILEDGEIGQILHDNVILKSIYYTGE
VNGTYYQFGKHYGNKKYRK

Host: Rabbit

Reactivity: Human

Applications: IF, WB-Ti, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 4675

Gene Symbol: NAP1L3

Gene Alias: MB20, MGC26312, NPL3

Gene Summary: This gene encodes a member of the nucleosome assembly protein (NAP) family. It is linked closely to a region of genes responsible for several X-linked mental retardation syndromes. Alternatively spliced transcript variants have been found for this gene, but the biological validity of the variants has not been determined. [provided by RefSeq]