

Datasheet

PDE1A MaxPab rabbit polyclonal antibody (D01)

Gene Symbol: PDE1A

Gene Alias: HCAM1, HSPDE1A, MGC26303

Catalog Number: H00005136-D01

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human PDE1A protein.

Immunogen: PDE1A (NP_005010.2, 1 a.a. ~ 545 a.a) full-length human protein.

Sequence:

MGSSATEIEELENTTFKYLTGEQTEKMWQRLKGILRCL
VKQLERGDVNVVDLKKNIEYAASVLEAVYIDETRRLD
TEDELSDIQTDSPSEVRDWLASTFTRKMGMTKKKPE
EKPKFRSIVHAVQAGIFVERMYRKTYHMGVGLAYPAAVI
VTLKDVDKWSFDVFALNEASGEHSLKFMIELFTRYDL
INRFKIPVSCLITFAEAELEVGYSKYKNPYHNLHAADV
QTVHYIMLHTGIMHWLTELEILAMVFAAAIHDEHTGTT
NNFHIQTRSDVAILYNDRSVLENHHVSAAYRLMQEEE
MNILINLSKDDWRDLRNLVIEMVLSTDMSGHFQQIKNIR
NSLQQPEGIDRAKTMSLILHAADISHPAKSWKLHYRWT
MALMEEFFLQGDKEAELGLPFSPKCDRKSTMVAQSQI
GFIDFIVEPTFSLTDSKIVPLIEEASKAETSSYVASS
STTIVGLHIADALRRSNTKGSMDSGSPDYSLAAVDL
KSFKNLVDIIQQNKERWKELAAQGESDLHKNSDLV
NAEEKHDETHS

Host: Rabbit

Reactivity: Human

Applications: IP, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 5136