

Datasheet

PGM3 MaxPab mouse polyclonal antibody (B01)

Gene Symbol: PGM3

Gene Alias: AGM1, DKFZp434B187, PAGM

Catalog Number: H00005238-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human PGM3 protein.

Immunogen: PGM3 (NP_056414.1, 1 a.a. ~ 542 a.a) full-length human protein.

Sequence:

MDLGAITKYSALHAKPONGLILQYGTAGFRTKAEHL DHV
MFRMGLLAVLR SKQTKSTIGVMVTASHNPEEDNGVKL
VDPLGEM LAPSWEEHATCLANAE EQDMQRVLIDISEK
EAVNLQQDAFVVIGRDTRPSSEKLSQSVIDGVTVLGG
QFHDYGLLTTPQLHYMVYCRNTGGRYGKATIEGYYQK
LSKAFVELTKQASCSGDEYRSLKVDCANGIGALKLRE
MEHYFSQGLSVQLFNDGSKGKLNHL CGADFKSHQK
PPQGMEIKSNERCCSFDGDADRIVYYH DADGHFHLI
DGDKIATLISSFLKELLVEIGESLNIGVVQTAYANGSSTR
YLEEVMKVPVYCTKTGVKHLHHKAQEFDIGVYFEANG
HGTALFSTAVEMKIQSAEQLEDKKRKA AKMLENIIDLF
NQAAGDAISDMLVIEAILALKGLTVQQWDALYTDLPNR
QLKVQVADRRVISTTDAERQAVTPPGLQE AINDLVKKY
KLSRAFVRPSGTEDVVRVYAEADSQESADHLAHEVSL
AVFQLAGGIGERPQPGF

Host: Mouse

Reactivity: Human

Applications: WB-Ti, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 5238