

Datasheet

POU2F1 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00005451-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human POU2F1 protein.

Immunogen: POU2F1 (AAH03571.1, 1 a.a. ~ 743 a.a) full-length human protein.

Sequence:

MNNPSETSKPSMESGDGNTGTQTNGLDQKQPVVVG
GAISTAQAQAFGLGHLHQVQLAGTSLQAAAQSLNVQSK
SNEESGDSQQPSQPSQQPSVQAAIPQTQLMLAGGQIT
GLTLTPAQQLLLQQAQAQQLLAAAVQQHSASQQH
SAAGATISASAATPMTQIPLSQPIQIAQDLQQLQQLQQ
QNLNLQQFVLVHPTTNLQPAQFIISQTPQGQQGLLQA
QNLLTQLPQQSQANLLQSQPSITLTSQPATPTRTIAAT
PIQTL PQSQSTPKRIDTPSLEEPSDLEEELEQFAKTFKQ
RRIKLGFTQGDVGLAMGKLYGNDFSQTTISRFEALNLS
FKNMCKLKPLLEKWLND AENLSSDSSLSSPSALNSPGI
EGLSRRRKRTSIETNIRVALEKSFLENQKPTSEEITMI
ADQLNMEKEVIRVWFCNRRQKEKRINPPSSGGTSSSP
IKAIFPSPTSLVATTPSLVTSSAATTLTVSPVLPLTSAAV
TNLSVTGTSDTTSNNTATVISTAPPASSAVTSPSLSPS
PSASASTSEASSASETSTTQTTSTPLSSPLGTSQVMVT
ASGLQTA AAAALQGAAQLPANASLAAMAAAAGLNPSL
MAPSQFAAGGALLSLNPGTL SGALSPALMSNSTLATIQ
ALASGGSLPITSLDATGNLVFANAGGAPNIVTAPLFLNP
QNLSLLTSNPVSLVSAAAASAGNSAPVASLHATSTSAE
SIQNSLFTVASASGAASSTTTASKAQ

Host: Mouse

Reactivity: Human

Applications: IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 5451

Gene Symbol: POU2F1

Gene Alias: OCT1, OTF1