

Datasheet

SAE1 polyclonal antibody

Catalog Number: PAB19246

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of SAE1.

Immunogen: A synthetic peptide corresponding to 20 amino acids near C-terminus of human SAE1.

Host: Rabbit

Reactivity: Human, Mouse, Rat

Applications: ELISA, WB-Ce

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Affinity purification

Concentration: 1 mg/mL

Recommend Usage: Western blot (1-2 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.02% sodium azide)

Storage Instruction: Store at 4°C for three months. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 10055

Gene Symbol: SAE1

Gene Alias: AOS1, FLJ3091, HSPC140, SUA1