

## Datasheet

### C1RL polyclonal antibody

**Gene Alias:** C1RL1, C1RLP, C1r-LP, CLSPa**Catalog Number:** PAB20548**Regulation Status:** For research use only (RUO)**Product Description:** Rabbit polyclonal antibody raised against recombinant C1RL.**Immunogen:** Recombinant protein corresponding to amino acids of human C1RL.**Sequence:**NVLPVCLPDNETLYRSGLLGYVSGFGMEMGWLTTELK  
YSRLPVAPREACNAWLQKRQRPEVFSDNMFCVGDET  
Q**Host:** Rabbit**Reactivity:** Human**Applications:** IHC-P, WB

(See our web site product page for detailed applications information)

**Protocols:** See our web site at<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols**Form:** Liquid**Purification:** Antigen affinity purification**Isotype:** IgG**Recommend Usage:** Western Blot (1:250-1:500)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)**Storage Instruction:** Store at 4°C. For long term storage store at -20°C.  
Aliquot to avoid repeated freezing and thawing.**Entrez GeneID:** 51279**Gene Symbol:** C1RL