

Datasheet

DOCK5 polyclonal antibody

Gene Symbol: DOCK5

Catalog Number: PAB20955

Gene Alias: DKFZp451J181, DKFZp779M164,
DKFZp781J211

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against recombinant DOCK5.

Immunogen: Recombinant protein corresponding to amino acids of human DOCK5.

Sequence:

RELKEKVEKHYGVITLPPNLTERKQSRTGSIVLPYIMSS
TLRRLSITSVTSSVVSTSSNSSDNAPSRPGSDGSILEPL
LERRASSGARVEDLSLREENSENRIKFKRKDWLSLK
SQVIAEKAPEPDLMSPTRKAQR

Host: Rabbit

Reactivity: Human

Applications: IF, IHC-P

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen affinity purification

Isotype: IgG

Recommend Usage: Immunohistochemistry

(1:200-1:500)

Immunofluorescence (1-4 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 80005