

Datasheet

CKAP5 polyclonal antibody

Gene Alias: CHTOG, FLJ35359, KIAA0097, MSPS, TOG, TOGp, ch-TOG

Catalog Number: PAB23651

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against recombinant CKAP5.

Immunogen: Recombinant protein corresponding to amino acids of human CKAP5.

Sequence:

MSSKLNQARSMGHPAAQMVRREFQLDLDEIENDN
GTVRCPELVLQHKLDDIFEPVLIPEPKIRAVSPHFDD
MHSNTASTINFI

Host: Rabbit

Reactivity: Human

Applications: IHC-P

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen affinity purification

Isotype: IgG

Recommend Usage: Immunohistochemistry

(1:1000-1:2500)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 9793

Gene Symbol: CKAP5