

Datasheet

UBXN6 polyclonal antibody

Gene Alias: DKFZp667D109, UBXD1, UBXDC2

Catalog Number: PAB23712

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against recombinant UBXN6.

Immunogen: Recombinant protein corresponding to amino acids of human UBXN6.

Sequence:

IYTFNKDQDRVKLGVDTI AKYLDNIHLHP EEEKYRKIKL
QNKVFQERINCLEGT HEFFEAIGFQKVLLPAQDQEDPE
EFYVLSET

Host: Rabbit

Reactivity: Human

Applications: IF, IHC-P

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen affinity purification

Isotype: IgG

Recommend Usage: Immunohistochemistry (1:20-1:50)

Immunofluorescence (1-4 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 80700

Gene Symbol: UBXN6