

## Datasheet

### MPV17L polyclonal antibody

**Catalog Number:** PAB23765

**Regulation Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against recombinant MPV17L.

**Immunogen:** Recombinant protein corresponding to amino acids of human MPV17L.

**Sequence:** SQQSGDGTFKSAFTILYTKGTSATEGYPKK

**Host:** Rabbit

**Reactivity:** Human

**Applications:** IHC-P

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Purification:** Antigen affinity purification

**Isotype:** IgG

**Recommend Usage:** Immunohistochemistry (1:10-1:20)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)

**Storage Instruction:** Store at 4°C. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 255027

**Gene Symbol:** MPV17L

**Gene Alias:** FLJ39599, MGC70356, MLPH1, MLPH2