



BSE / TSE Statement

VWR* TRI-POUR BEAKERS

This document certifies that the product identified below meets the following criteria:

North American Catalog N°:	25384-156
Description:	VWR BEAKER GRAD DSP 50ML CS100

The product may contain North American cattle derived additives. Any tallow used has been treated by transesterification, or hydrolysis at 200 degrees C for 20 minutes at 40 bar, or greater.

Signed:

A handwritten signature in blue ink, appearing to read "Ken Crossley", with a stylized flourish at the end.

Ken Crossley
Manager
Quality Assurance

Date: July 27, 2011