

Datasheet

TOLLIP polyclonal antibody

Catalog Number: PAB9721

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of TOLLIP.

Immunogen: A synthetic peptide corresponding to C-terminus of human TOLLIP.

Host: Rabbit

Reactivity: Human, Mouse, Rat

Applications: IHC, WB-Ti

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.02% sodium azide)

Storage Instruction: Store at -20°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 54472

Gene Symbol: TOLLIP

Gene Alias: FLJ33531, IL-1RAcPI

References:

1. The expanding family of MyD88-like adaptors in Toll-like receptor signal transduction. McGettrick AF, O'Neill LA. Mol Immunol. 2004 Jul;41(6-7):577-82.
2. Toll-like receptors. Takeda K, Kaisho T, Akira S. Annu Rev Immunol. 2003;21:335-76. Epub 2001 Dec 19.