

## Datasheet

### TOLLIP polyclonal antibody

**Catalog Number:** PAB9722

**Regulation Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of TOLLIP.

**Immunogen:** A synthetic peptide corresponding to internal region of human TOLLIP.

**Host:** Rabbit

**Reactivity:** Human, Mouse, Rat

**Applications:** IHC, WB-Ti

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Recommend Usage:** The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS (0.02% sodium azide)

**Storage Instruction:** Store at -20°C.  
Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 54472

**Gene Symbol:** TOLLIP

**Gene Alias:** FLJ33531, IL-1RAcPIP