

Datasheet

C10orf116 polyclonal antibody

Catalog Number: PAB9766

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of C10orf116.

Immunogen: A synthetic peptide corresponding to human C10orf116.

Host: Rabbit

Reactivity: Human

Applications: WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.08% sodium azide)

Storage Instruction: Store at -20°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 10974

Gene Symbol: C10orf116

Gene Alias: APM2

Gene Summary: APM2 gene is exclusively expressed in adipose tissue. Its function is currently unknown.
[provided by RefSeq]