

Datasheet

BCL10 polyclonal antibody

Catalog Number: PAB9776

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of BCL10.

Immunogen: A synthetic peptide corresponding to amino acids 5-19 of human BCL10.

Host: Rabbit

Reactivity: Human

Applications: WB-Ce

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.02% sodium azide)

Storage Instruction: Store at -20°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 8915

Gene Symbol: BCL10

Gene Alias: CARMEN, CIPER, CLAP, c-E10, mE10

References:

1. Equine herpesvirus-2 E10 gene product, but not its cellular homologue, activates NF-kappaB transcription factor and c-Jun N-terminal kinase. Thome M, Martinon F, Hofmann K, Rubio V, Steiner V, Schneider P, Mattmann C, Tschopp J. J Biol Chem. 1999 Apr 9;274(15):9962-8.
2. CIPER, a novel NF kappaB-activating protein containing a caspase recruitment domain with homology to Herpesvirus-2 protein E10. Koseki T, Inohara N, Chen

S, Carrio R, Merino J, Hottiger MO, Nabel GJ, Nunez G. J Biol Chem. 1999 Apr 9;274(15):9955-61.

3. Bcl10 is involved in t(1;14)(p22;q32) of MALT B cell lymphoma and mutated in multiple tumor types. Willis TG, Jadayel DM, Du MQ, Peng H, Perry AR, Abdul-Rauf M, Price H, Karran L, Majekodunmi O, Wlodarska I, Pan L, Crook T, Hamoudi R, Isaacson PG, Dyer MJ. Cell. 1999 Jan 8;96(1):35-45.