

Datasheet

Bub1b polyclonal antibody

Catalog Number: PAB9783

Regulation Status: For research use only (RUO)

Product Description: Sheep polyclonal antibody raised against partial recombinant Bub1b.

Immunogen: Recombinant protein corresponding to amino acids 232-603 of mouse Bub1b.

Host: Sheep

Reactivity: Mouse

Applications: IHC, IP

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.08% sodium azide)

Storage Instruction: Store at -20°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 12236

Gene Symbol: Bub1b

Gene Alias: AU045529, BUBR1, MGC36092

References:

1. Mutation and expression analysis of human BUB1 and BUB1B in aneuploid breast cancer cell lines. Myrie KA, Percy MJ, Azim JN, Neeley CK, Petty EM. Cancer Lett. 2000 May 1;152(2):193-9.
2. The mouse mitotic checkpoint gene bub1b, a novel bub1 family member, is expressed in a cell cycle-dependent manner. Davenport JW, Fernandes ER, Harris LD, Neale GA, Goorha R. Genomics. 1999 Jan 1;55(1):113-7.