

Datasheet

Tlr7 polyclonal antibody

Catalog Number: PAB9788

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of Tlr7.

Immunogen: A synthetic peptide corresponding to internal region of mouse Tlr7.

Host: Rabbit

Reactivity: Human, Mouse

Applications: ICC, WB-Ce

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.02% sodium azide)

Storage Instruction: Store at -20°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 170743

Gene Symbol: Tlr7

Gene Alias: -

References:

1. The expanding family of MyD88-like adaptors in Toll-like receptor signal transduction. McGettrick AF, O'Neill LA. Mol Immunol. 2004 Jul;41(6-7):577-82.
2. Ligand-regulated chimeric receptor approach reveals distinctive subcellular localization and signaling properties of the Toll-like receptors. Nishiya T, DeFranco AL. J Biol Chem. 2004 Apr 30;279(18):19008-17. Epub 2004 Feb 19.