

Datasheet

CHEK1 polyclonal antibody

Catalog Number: PAB9841

Regulation Status: For research use only (RUO)

Product Description: Sheep polyclonal antibody raised against partial recombinant CHEK1.

Immunogen: Recombinant GST fusion protein corresponding to amino acids 1-474 of human CHEK1.

Host: Sheep

Reactivity: Human, Monkey, Mouse, Rat

Applications: IP, WB-Ce

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Recommend Usage: Western Blot (1-10 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.08% sodium azide)

Storage Instruction: Store at -20°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 1111

Gene Symbol: CHEK1

Gene Alias: CHK1