

## Datasheet

### Fli1 polyclonal antibody

**Catalog Number:** PAB9876

**Regulation Status:** For research use only (RUO)

**Product Description:** Sheep polyclonal antibody raised against partial recombinant Fli1.

**Immunogen:** Recombinant GST fusion protein corresponding to amino acids 1-239 of mouse Fli1.

**Host:** Sheep

**Reactivity:** Human, Mouse

**Applications:** WB-Ce

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Recommend Usage:** Western Blot (1:500)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS (0.08% sodium azide)

**Storage Instruction:** Store at -20°C.

Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 14247

**Gene Symbol:** Fli1

**Gene Alias:** EWSR2, Fli-1, SIC-1, Sic1