

Datasheet

SYCP1 polyclonal antibody

Catalog Number: PAB12015

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of SYCP1.

Immunogen: A synthetic peptide corresponding to C-terminus of human SYCP1.

Host: Rabbit

Reactivity: Human, Mouse, Rat

Applications: IF, IHC-P

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody is specific to human SCP1.

Form: Liquid

Recommend Usage: Immunofluorescence

(1:300-1:750)

The optimal working dilution should be determined by the end user.

Storage Buffer: In Tris-citrate/phosphate buffer, pH 7.0-8.0 (0.1% sodium azide)

Storage Instruction: Store at 4°C. Do not freeze.

Entrez GeneID: 6847

Gene Symbol: SYCP1

Gene Alias: HOM-TES-14, MGC104417, SCP1

References:

1. Heritable and stable gene knockdown in rats. Dann CT, Alvarado AL, Hammer RE, Garbers DL. Proc Natl Acad Sci U S A. 2006 Jul 25;103(30):11246-51. Epub 2006 Jul 14.