

Datasheet

Sstr1 polyclonal antibody

Catalog Number: PAB12019

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of Sstr1.

Immunogen: A synthetic peptide corresponding to C-terminus of rat Sstr1.

Host: Rabbit

Reactivity: Human, Mouse, Rat

Applications: ICC, IP, WB-Ti

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Recommend Usage: Western Blot (1:1000-1:5000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In antiserum (0.05% sodium azide)

Storage Instruction: Store at -20°C or -80°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 25033

Gene Symbol: Sstr1

Gene Alias: Gpcrrna

References:

1. Development and use of a receptor antibody to characterize the interaction between somatostatin receptor subtype 1 and G proteins. Gu YZ, Brown PJ, Loose-Mitchell DS, Stork PJ, Schonbrunn A. Mol Pharmacol. 1995 Dec;48(6):1004-14.