

Datasheet

Aggf1 polyclonal antibody

Catalog Number: PAB12151

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of Aggf1.

Immunogen: A synthetic peptide corresponding to amino acids 300-400 of mouse Aggf1.

Host: Rabbit

Reactivity: Mouse

Applications: IHC, WB-Ti

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody is specific to mouse VG5Q.

Form: Liquid

Recommend Usage: Immunohistochemistry (1:250)

Western Blot (1:500)

The optimal working dilution should be determined by the end user.

Storage Buffer: In antiserum

Storage Instruction: Store at 4°C for short term. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 66549

Gene Symbol: Aggf1

Gene Alias: 2010009L17Rik, 2310029P06Rik, AW112072, Peg3, VG5Q

References:

1. Identification of an angiogenic factor that when mutated causes susceptibility to Klippel-Trenaunay syndrome. Tian XL, Kadaba R, You SA, Liu M, Timur

AA, Yang L, Chen Q, Szafranski P, Rao S, Wu L, Housman DE, DiCorleto PE, Driscoll DJ, Borrow J, Wang Q. Nature. 2004 Feb 12;427(6975):640-5.