

Datasheet

Adfp polyclonal antibody

Catalog Number: PAB12152

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of Adfp.

Immunogen: A synthetic peptide corresponding to amino acids 350-425 of mouse Adfp.

Host: Rabbit

Reactivity: Human, Mouse

Applications: ICC, IF, WB-Ti

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Recommend Usage: Immunocytochemistry

(1:50-1:200)

Immunofluorescence (1:50-1:200)

Western Blot (2 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (30% glycerol)

Storage Instruction: Store at -20°C or -80°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 11520

Gene Symbol: Adfp

Gene Alias: AA407157, ADPH, Adrp

References:

1. Fixation methods for the study of lipid droplets by immunofluorescence microscopy. DiDonato D, Brasaemle DL. J Histochem Cytochem. 2003 Jun;51(6):773-80.