

Datasheet

REPIN1 polyclonal antibody (A01)

Catalog Number: H00029803-A01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a partial recombinant REPIN1.

Immunogen: REPIN1 (NP_037532, 1 a.a. ~ 110 a.a)
partial recombinant protein with GST tag.

Sequence:

MLERRCRGPLAMGLAQPRLLSGPSQESPQTLGKESR
GLRQQGTSVAQSGAQAPGRAHRCACHRRHFPGWVA
LWLHTRRCQARLPLPCPEGRRFRHAPFLALHRQVH
AAA

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 29803

Gene Symbol: REPIN1

Gene Alias: AP4, RIP60, ZNF464, Zfp464