

Datasheet

EHD3 monoclonal antibody (M01), clone 4B7

Catalog Number: H00030845-M01

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant EHD3.

Clone Name: 4B7

Immunogen: EHD3 (NP_055415, 357 a.a. ~ 406 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

KRMQDQLQAQDFSKFQPLKSKLLEVDDMLAHDIAQL
MVLVRQEESQRPI

Host: Mouse

Reactivity: Human

Applications: ELISA, IF, IHC-P, S-ELISA, WB-Ce,
WB-Re

(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG1 Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 30845

Gene Symbol: EHD3

Gene Alias: PAST3