

Datasheet

TAX1BP3 monoclonal antibody (M01), clone 4A10

Borthwick K, Donne AJ, Hampson L, Hampson IN. Br J Cancer. 2010 Dec 7. [Epub ahead of print]

Catalog Number: H00030851-M01

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a full length recombinant TAX1BP3.

Clone Name: 4A10

Immunogen: TAX1BP3 (AAH23980, 1 a.a. ~ 124 a.a)
full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

MSYIPGQPVTAVVQRVEIHKLRQGENLILGFSIGGGIDQ
DPSQNPFSKTDKGIYVTRVSEGGPAEIALQLIGDKI
MQVNGWDMTMVTHDQARKRLTKRSEEVVRLLVTRQS
LQKAVQQSMLS

Host: Mouse

Reactivity: Human

Applications: ELISA, S-ELISA, WB-Re
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG1 kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 30851

Gene Symbol: TAX1BP3

Gene Alias: TIP-1

References:

1. The HPV16 E6 binding protein Tip-1 interacts with ARHGEF16, which activates Cdc42. Oliver AW, He X,