

Datasheet

GMNN monoclonal antibody (M02), clone 1B10

Catalog Number: H00051053-M02

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant GMNN.

Clone Name: 1B10

Immunogen: GMNN (AAH05185, 110 a.a. ~ 209 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

LYEALKENEKLHKEIEQKDNEIARLKKENKELAEVAEHV
QYMAELIERLNGEPLDNFESLDNQEFDSEEETVEDSLV
EDSEIGTCAEGTVSSSTDAKPCI

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re, WB-Tr
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 51053

Gene Symbol: GMNN

Gene Alias: Gem, RP3-369A17.3