

Datasheet

LAP3 polyclonal antibody (A01)

Catalog Number: H00051056-A01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a partial recombinant LAP3.

Immunogen: LAP3 (NP_056991, 420 a.a. ~ 519 a.a)
partial recombinant protein with GST tag.

Sequence:

EASLETGDRVWRMPLFEHYTRQVVDCQLADVNNIGKY
RSAGACTAAAFLEFVTHPKWAHLDIAGVMTNKDEVP
YLRKGMTGRPTRTLIEFLRFSQDNA

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Ce, WB-Re, WB-Tr
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 51056

Gene Symbol: LAP3

Gene Alias: LAP, LAPEP, PEPS

References:

1. Cytosolic Aminopeptidases Influence MHC Class
I-Mediated Antigen Presentation in an Allele-Dependent
Manner. Kim E, Kwak H, Ahn K. J Immunol. 2009 Dec
1;183(11):7379-87. Epub 2009 Nov 16.