

Datasheet

SH3GLB1 monoclonal antibody (M01), clone 1B3-A5

Gene Symbol: SH3GLB1

Gene Alias: Bif-1, CGI-61, KIAA0491, dJ612B15.2

Catalog Number: H00051100-M01

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a full length recombinant SH3GLB1.

Clone Name: 1B3-A5

Immunogen: SH3GLB1 (AAH07455, 1 a.a. ~ 365 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

MNIMDFNVKKLAADAGTFLSRAVQFTEEKLGQAEKTE
LDAHLENLLSKAECTKIWTEKIMKQTEVLLQPNPNARIE
EFVYEKLDRKAPSRINNPELLGQYMIDAGTEFGPGTAY
GNALIKCGETQKRIGTADRELIQTSALNFLTPLRNFIEG
DYKTIAKERKLLQNKRLDLDAKTRLKKAKAAETRNS
EQELRITQSEFDRQAEITRLLLEGISSTHAHHLRCLNDF
VEAQMTYYAQCYQYMLDLQKQLGSFSPNYLSNNNQ
SVTPVPSVLPNAIGSSAMASTSGLVITSPSNLSDLKEC
SGSRKARVLYDYDAANSTELSLLADEVITVFSVVGMD
DWLMGERGNQKGKVPITYLELLN

Host: Mouse

Reactivity: Human

Applications: ELISA, IHC-P, S-ELISA, WB-Ce, WB-Re, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG1 kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 51100