

Datasheet

NDUFAF1 MaxPab mouse polyclonal antibody (B01)

OMIM]

Catalog Number: H00051103-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human NDUFAF1 protein.

Immunogen: NDUFAF1 (NP_057097, 1 a.a. ~ 327 a.a) full-length human protein.

Sequence:

MALVHKLLRGTYFLRKFSKPTSALYPFLGIRFAEYSSSL
QKPVASPGKASSQRKTEGDLQGDHQKEVALDITSSEE
KPDVSFDKAIRDEAIYHFRLKDEIVDHWRGPEGHPLH
EVLLEQAKVWVWQFRGKEDLDKWTVTSDKTIGGRSEVF
LKMKGKNNQSALLYGTLSSSEAPQDGESTRSGYCAMISR
IPRGAFERKMSYDWSQFNTLYLRVRGDGRPWMVNIK
EDTDFQRTNQMYSYFMFTRGGPYWQEVKIPFSKFFF
SNRGRIRDVQHELPLDKISSIGFTLADKVDGPFFLEIDFI
GVFTDPAHTEEFAYENSPELNPRLFK

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 51103

Gene Symbol: NDUFAF1

Gene Alias: CGI-65, CGI65, CIA30

Gene Summary: The NDUFAF1 gene encodes the human homolog of a *Neurospora crassa* Cia30 protein involved in the assembly of complex I.[supplied by