

Datasheet

MRT04 MaxPab rabbit polyclonal antibody (D01)

Catalog Number: H00051154-D01

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human MRT04 protein.

Immunogen: MRT04 (NP_057267.2, 1 a.a. ~ 239 a.a) full-length human protein.

Sequence:

MPKSKRDKKVSLTKTAKKGLELKQNLIEELRKCVDTYK
YLFIFSVANMRNSKLKDIRNAWKHSRMFFGKNKVMV
ALGRSPSDEYKDNLHQVSKRLRGEVGLLFTNRTKEEV
NEWFTKYTEM DYARAGNKA AFTVSLDPGPLEQFPHS
MEPQLRQLGLPTALKRGVVTLSDYEVCKEGDVLTP
QARVLKLFGEYMAEFKVTIKYMWDSQSGRFQMGDD
LPESASESTEESDSEDD

Host: Rabbit

Reactivity: Human

Applications: IP, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 51154

Gene Symbol: MRT04

Gene Alias: C1orf33, MRT4, dJ657E11.4

Gene Summary: This gene encodes a protein sharing a low level of sequence similarity with ribosomal protein P0. While the precise function of the encoded protein is currently unknown, it appears to be involved in mRNA turnover and ribosome assembly. [provided by RefSeq]