

## Datasheet

### SHISA5 monoclonal antibody (M03), clone 3G5

**Catalog Number:** H00051246-M03

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody  
raised against a full-length recombinant SHISA5.

**Clone Name:** 3G5

**Immunogen:** SHISA5 (AAH01463, 1 a.a. ~ 137 a.a)  
full-length recombinant protein with GST tag. MW of the  
GST tag alone is 26 KDa.

**Sequence:**

MGFGATLAVGLTIFVLSVVTIIICFTSCCCLYKTCRRP  
RPVVTTTTSTTVVHAPYPQPPSVPPSYGPSYQGYHT  
MPPQPGMPAAPYPMPQYPPPYPAQPMGPPAYHETLA  
GGAAAPYPASQPPYNPAYMDAPKAAL

**Host:** Mouse

**Reactivity:** Human

**Applications:** ELISA, WB-Re

(See our web site product page for detailed applications  
information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product  
page for detailed protocols

**Isotype:** IgG2a Kappa

**Storage Buffer:** In 1x PBS, pH 7.4

**Storage Instruction:** Store at -20°C or lower. Aliquot to  
avoid repeated freezing and thawing.

**Entrez GeneID:** 51246

**Gene Symbol:** SHISA5

**Gene Alias:** SCOTIN, hShisa5