

Datasheet

PPP5C MaxPab rabbit polyclonal antibody (D01)

Gene Alias: FLJ36922, PP5, PPP5

Catalog Number: H00005536-D01

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human PPP5C protein.

Immunogen: PPP5C (NP_006238.1, 1 a.a. ~ 499 a.a) full-length human protein.

Sequence:

MAMAEGERTECAEPPRDEPPADGALKRAEELKTQAN
DYFKAQDYENAIKFYSQAIELNPSNAIYYGNRSLAYLRT
ECYGYALGDATRAIELDKKYIKGYRRAASNMALGKFR
AALRDYETVVKVPHDKDAKMKYQECNKIVKQKAER
AIAGDEHKRSVVDSLDESMTIEDEYSGPKLEDGKVTIS
FMKELMQWYKDQKKLHRKCAYQILVQVKEVLSKLSTL
VETTLKETEKITVCGDTHGQFYDLLNIFELNGLPSETNP
YIFNGDFVDRGSFSVEVILTFLGFKLLYPDFHLLRGNH
ETDNMNQIYGFEGEVKAKYTAQMYELFSEVFEWLPLA
QCINGKVLIMHGGLFSEDGVTLDIRKIERNRQPPDSG
PMCDLLWSDPQPQNGRSISKRGVSCQFGPDVTKAFL
EENLDYIIRSHEVKAEGYEVAHGGRCVTVFSAAPNYC
DQMGNKASYIHLQGSDLRPQFHQFTAVPHPNVKPMA
YANTLLQLGMM

Host: Rabbit

Reactivity: Human

Applications: IP, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 5536

Gene Symbol: PPP5C