

Datasheet

PRL monoclonal antibody (M06), clone X1

Catalog Number: H00005617-M06

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a full length recombinant PRL.

Clone Name: X1

Immunogen: PRL (AAH15850, 29 a.a. ~ 227 a.a)
full-length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

LPICPGGAARCQVTLRDLFDRAVVLSHYIHNLSEMFS
EFDKRYTHGRGFITKAINCHTSSLATPEDKEQAQQM
NQKDFLSLIVSILRSWNEPLYHLVTEVRGMQEAPAILS
KAVEIEEQTKRLLEGMEIVSQVHPETKENEIYPVWSG
LPSLQMADEESRLSAYYNLLHCLRRDSHKIDNYLKKL
CRIIHNNNC

Host: Mouse

Reactivity: Human

Applications: ELISA, S-ELISA, WB-Re, WB-Tr
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 5617

Gene Symbol: PRL

Gene Alias: -