

Rabbit Polyclonal Antibody to Human CD117/c-kit

Catalog No.:	PDR 045
Applications:	Analyte Specific Reagent. Analytical and performance characteristics are not established. Optimal working conditions must be determined by the end user.
Immunogen:	A synthetic peptide from c-terminus of human CD117/c-kit protein.
Host:	Rabbit
Specificity:	This antibody reacts with a 145 kD protein. c-kit is a transmembrane receptor tyrosine kinase that is expressed in many tissues and cells. c-kits is involved in the development of several lineages of stem cells, such as germ cells, neural crest derived melanocytes, and hematopoietic precursor cells.
Positive Control:	Tonsil
Cellular Localization:	Cell membrane, cytoplasm
Storage:	Store at 2-8°C. Do not use beyond the expiration date stated on the label.
References:	i) Arber et al. Hum Pathol 29: 498, 1998. ii) Rajpert et al. J Pathol 178: 166, 1996.

ASR: Analyte Specific Reagent

DBS will not be held responsible for patent infringement or other violation that may occur with the use of our product

DBS

1020 Serpentine Lane, # 114, Pleasanton, CA 94566 Tel: 925 484 3350, Fax: 925 484 3390

Website: www.dbiosys.com e-mail: customersupport@dbiosys.com