

Datasheet

SKIL purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00006498-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human SKIL protein.

Immunogen: SKIL (NP_005405.1, 1 a.a. ~ 684 a.a) full-length human protein.

Sequence:

MENLQTNFSLVQGSTKKLNGMGDDGSPPAKKMITDIH
VNGKTINKVPTVKKEHLDDYGEAPVETDGEHVKRTCT
SVPETLHLNPSLKHTLAQFHLSSQSSLGGPAAFSARH
SQESMSPTVFLPLSPQVLPGLLIPSDSSTELTQTVL
EGESISCFQVGGEKRLCLPQVLNSVLREFTLQQINTVC
DELYIYCSRCTSDQLHILKVLGILPFNAPSCGLITLTDQA
RLCNALLRPRTFPQNGSVLPAKSSLAQLKETGSAFEV
EHECLGKCQGLFAPQFYVQPDAPCIQCCECCGMFAP
QTFVMHSHRSPDKRTCHWGFESAKWHCYLHVNQKY
LGTPEEKKLKIILEEMKEKFSMRSGKRNQSKTDAPSG
MELQSWYPVIKQEGDHVSQTHSFLHPSYYLYMCDKV
VAPNVSLTSAVSQSKELTKTEASKSISRQSEKAHSSGK
LQKTVSYPDVSLEEQEKMDLKTSRELCSRLDASISNNS
TSKRKSESATCNLVRDINKVGIGLVAAASSPLLVKDVIC
EDDKGKIMEEVMRTYLKQKEKLNILQKKQQLQMEVK
MLSSSKSMKELTEEQNLQKELESLQNEHAQRMEEF
YVEQKDLEKKLEQIMKQKCTCDSNLEKDKEAEYAGQL
AELRQRDLHAEADRQELQDELQREARQKLEMMIKE
LKLQILKSSKTAKE

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to

avoid repeated freezing and thawing.

Entrez GeneID: 6498

Gene Symbol: SKIL

Gene Alias: SNO, SnoA, SnoN