

Datasheet

TNNT3 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00007140-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human TNNT3 protein.

Immunogen: TNNT3 (NP_006748.1, 1 a.a. ~ 250 a.a) full-length human protein.

Sequence:

MSDEEVEQVEEQYEEEEEAQEEEEVQEEKPRPKLT
APKIPEGEKVDFDDIQKKRQNKDLMELQALIDSHFEAR
KKEEEELVALKERIEKRRAERAEQQRIRAEKERERQN
RLAEEKARREEEDAKRRAEDDLKKKKALSSMGANYSS
YLAKADQKRGKKQTAREMKKKILAERRKPLNIDHLGED
KLRDKAKELWETLHQLEIDKFEFGEKLRQKYDITTLR
SRIDQAQKHSKKAGTPAKGKVGGRWK

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 7140

Gene Symbol: TNNT3

Gene Alias: AMCD2B, DA2B, DKFZp779M2348, FSSV