

Datasheet

TRPC1 monoclonal antibody (M03), clone 1E4

Catalog Number: H00007220-M03

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant TRPC1.

Clone Name: 1E4

Immunogen: TRPC1 (NP_003295, 442 a.a. ~ 505 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

AHNKFHDFADRKDWDAFHPTLVAEGLFAFANVLSYLR
LFFMYTTSSILGPLQISMGMQLQDFGK

Host: Mouse

Reactivity: Human

Applications: ELISA, S-ELISA

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2b Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 7220

Gene Symbol: TRPC1

Gene Alias: HTRP-1, MGC133334, MGC133335, TRP1