

## Datasheet

### F2RL3 polyclonal antibody

**Catalog Number:** PAB26183

**Regulation Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of F2RL3.

**Immunogen:** A synthetic peptide corresponding to 17 amino acid at cytoplasmic domain of human F2RL3.

**Host:** Rabbit

**Reactivity:** Human

**Applications:** IHC-P

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** BLAST analysis of the peptide immunogen showed no homology with other human proteins.

**Form:** Liquid

**Purification:** Immunoaffinity chromatography

**Recommend Usage:** Immunohistochemistry

(Formalin/PFA-fixed paraffin-embedded sections) (5-10 ug/mL)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS (0.1% sodium azide)

**Storage Instruction:** Store at 4°C. For long term storage store at -80°C.

Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 9002

**Gene Symbol:** F2RL3

**Gene Alias:** PAR4

**Gene Summary:** Coagulation factor II (thrombin)

receptor-like 3 (F2RL3) is a member of the large family of 7-transmembrane-region receptors that couple to guanosine-nucleotide-binding proteins. F2RL3 is also a member of the protease-activated receptor family. F2RL3 is activated by proteolytic cleavage of its extracellular amino terminus. The new amino terminus functions as a tethered ligand and activates the receptor. F2RL3 is activated by thrombin and trypsin. [provided by RefSeq]