

Datasheet

GPR18 polyclonal antibody

Catalog Number: PAB26338

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of GPR18.

Immunogen: A synthetic peptide corresponding to 17 amino acids at 2nd extracellular domain of human GPR18.

Host: Rabbit

Reactivity: Gorilla, Human, Monkey

Applications: IHC-P

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: BLAST analysis of the peptide immunogen showed no homology with other human proteins.

Form: Liquid

Purification: Immunoaffinity chromatography

Recommend Usage: Immunohistochemistry

(Formalin/PFA-fixed paraffin-embedded sections) (2.5 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.1% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -80°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 2841

Gene Symbol: GPR18

Gene Alias: -