

PDS No. 3xx311	PRODUCT DATA SHEET				Page 1 of 1
Revision 02	PE Disc Cap				 greiner bio-one
	Greiner Item-No. 3xx311				
Valid for Item-No.:	301311	302311	303311	322311	

1.	Description / Specification	
1.1	Description	PE disc cap 301311: for tube Art.-No. 103xxx 302311: for tube Art.-No. 106xxx, 109xxx 303311: for tube Art.-No. 112xxx, 115xxx 322311: for tube Art.-No. 199xxx
1.2	Dimensions	301311: External diameter: 10,8 mm 302311: External diameter: 12,1 mm 303311: External diameter: 13,3 mm 322311: External diameter: 24,0 mm
1.3	Volume	N/A
1.4	Material / Resin	LDPE (Low Density Polyethylene), free of heavy metal
1.5	Colour	Translucent
1.6	Sterilization	No
1.7	Quality Control	<u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	-
2.2	Temperature range	-50°C to +80°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: www.gbo.com/bioscience → Products → Literature → Technical Information → Chemical Resistance of Resins
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	301311	302311	303311	322311
3.1	Pieces / Bag	1000	500	500	7500
3.2	Pieces / Box	35000	30000	20000	7500
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)			
3.4	Other Information	-			

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 01	Date 7 January 2011	Date 24 January 2011	Date 25 January 2011	
Date 11.06.2008	Name S. Kaelberer	Name Dr. A. Fuhrmann	Name A. Schulz	