

PDS No. 6732xx	PRODUCT DATA SHEET					Page 1 of 1
Revision 03	8-tube PCR Strips, 0,2 ml, without Cap					 greiner bio-one
	Greiner Item-No. 6732xx					
Valid for Item-No.:	673210	673271	673272	673273	673274	
	673275	673276	673277			

1.	Description / Specification	
1.1	Description	8-tube PCR strips, 0,2 ml, without cap
1.2	Dimensions	<u>Tube</u> : length: 20,83 mm outside diameter: 7,09 mm inner diameter: 5,36 mm <u>Tube strip</u> : length: 71,12 mm
1.3	Volume	<u>Tube</u> : 0,2 ml
1.4	Material / Resin	PP (Polypropylene), free of heavy metal
1.5	Colour	673210: translucent 673271: assorted (orange, red, blue, green, yellow, violet) 673272: orange 673273: red 673274: blue 673275: green 673276: yellow 673277: violet
1.6	Sterilization	No
1.7	Quality Control	<u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	- For single use only - Suitable for Item-No. 37327x / 373281 and 373250

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens
2.2	Temperature range	-80°C to +105°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A (with a special centrifuge holder)
2.5	Chemical Resistance	See homepage: www.gbo.com/bioscience →Products →Literature →Technical Information →Chemical Resistance of Resins
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	673210, -272, -273, -274, -275, -276, -277	673271 (assorted)
3.1	Pieces / Bag	125	126
3.2	Pieces / Box	1250	1260
3.3	Lot-No.	YYMMDDNNN (year, month, day, mould number)	YYMMDDNNN (year, month, day, mould number)
3.4	Other Information	-	-

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 02	Date 12 March 2014	Date 13 March 2014	Date 13 March 2014	
Date 18.12.2009	Name S. Kaelberer	Name Dr. R. Heller	Name A. Schulz	