

PDS No. 6832xx	PRODUCT DATA SHEET					Page 1 of 1
Revision 04	Thin Wall PCR Tubes, 0,2 ml, Flat Cap					 greiner bio-one
	Greiner Item-No. 6832xx					
Valid for Item-No.:	683201	683271	683272	683273	683274	
	683275	683276	683277			

1.	Description / Specification	
1.1	Description	PP thin wall PCR tubes, 0,2 ml with attached flat cap
1.2	Dimensions	Length (without cap): 20,75 mm Outside diameter (tube): 7,09 mm Inner diameter (tube): 5,38 mm
1.3	Volume	0,2 ml
1.4	Material / Resin	PP (Polypropylene), free of heavy metal
1.5	Colour	683201: translucent 683271: assorted (orange, red, blue, green, yellow, violet) 683272: orange 683273: red 683274: blue 683275: green 683276: yellow 683277: violet
1.6	Sterilization	No
1.7	Quality Control	<u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens
2.2	Temperature range	-80°C to +105°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	4000 x g: fixed-angle rotor
2.5	Chemical Resistance	See homepage: www.gbo.com/bioscience →Products →Literature →Technical Information →Chemical Resistance of Resins
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	500
3.2	Pieces / Box	1000
3.3	Lot-No.	YYMMDDNNN (year, month, day, mould number)
3.4	Other Information	-

4.	Other Information	
	-	

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 03	Date 17 December 2009	Date 18 December 2009	Date 18 December 2009	
Date 16.11.2005	Name S. Kaelberer	Name Dr. L. Marchetti	Name A. Schulz	