

PDS No. 684201	PRODUCT DATA SHEET	Page 1 of 1
Revision 02	Thin Wall PCR Tubes, 0,2 ml, without Cap	 greiner bio-one
	Greiner Item-No. 684201	

1.	Description / Specification	
1.1	Description	Thin wall PCR tubes, 0,2 ml without cap
1.2	Dimensions	Length: 20,83 mm Outside diameter: 7,10 mm Inner diameter: 5,36 mm
1.3	Volume	0,2 ml
1.4	Material / Resin	PP (Polypropylene), free of heavy metal
1.5	Colour	Translucent
1.6	Sterilization	No
1.7	Quality Control	<u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens
2.2	Temperature range	-80°C to +105°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	4000 x g: fixed-angle rotor
2.5	Chemical Resistance	See homepage: www.gbo.com/bioscience →Products →Literature →Technical Information →Chemical Resistance of Resins
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	500
3.2	Pieces / Box	1000
3.3	Lot-No.	YYMMDDNNN (year, month, day, mould number)
3.4	Other Information	-

4.	Other Information	
	-	

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 01	Date 17 December 2009	Date 18 December 2009	Date 18 December 2009	
Date 06.12.2005	Name S. Kaelberer	Name Dr. R. Heller	Name A. Schulz	