

# SAFETY DATA SHEET

Cat# K163-100, CaspSELECT™ Caspase-3 Immunoassay Kit

SDS DATE: 13APR2015

## SECTION 1: PRODUCT AND COMPANY IDENTIFICATION

**PRODUCT NAME:** CaspSELECT™ Caspase-3 Immunoassay Kit  
**PRODUCT CODES:** Cat# K163-100  
**MANUFACTURER:** BioVision, Inc.  
**DIVISION:**  
**ADDRESS:** 155 S. Milpitas Blvd. Milpitas, CA 95035  
**EMERGENCY PHONE:** 858-373-8066  
**CHEMTREC PHONE:**  
**OTHER CALLS:** 408-493-1800  
**FAX PHONE:** 408-493-1801

## SECTION 2: HAZARDS IDENTIFICATION

Component	Description	Volume	Safety Information
Cell Lysis Buffer	Proprietary Buffer	25 ml	See below
Coating Buffer	Proprietary Buffer	10 ml	No hazards
Anti-Caspase-3 Antibody (20X)	liquid	0.5 ml	No hazards
Blocking Buffer	Proprietary Buffer	15 ml	No hazards
Incubation Buffer	Proprietary Buffer	100 ml	See below
DTT	Solution in H <sub>2</sub> O (1 M)	400 µl	See below
DEVD-AFC Substrate	In DMSO	500 µl	See below
Positive Control	Caspase-3, human recombinant	10 units	No hazards
Microtiter Plate	--	12 x 8 wells	No hazards
Adhesive Plate Cover	--	2 sheets	No hazards

### Triton X-100:

#### Emergency Overview

**OSHA Hazards:** Irritant, Harmful by ingestion

**GHS Classification:**  
Serious eye damage (Category 1)  
Skin irritation (Category 3)  
Acute aquatic toxicity (Category 3)  
Chronic aquatic toxicity (Category 3)

#### GHS Label elements, including precautionary statements

**Pictogram:**



**Signal word:**

Warning

**Hazard statement(s):**

H302: Harmful if swallowed  
H316 Causes mild skin irritation.  
H320: Causes eye irritation

**Precautionary statement(s):** P261 Avoid breathing dust/fumes/gas/mist/vapors/spray.

P264 Wash hands thoroughly after handling.

P273 Avoid release to the environment.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P403+P233 Store in a well-ventilated place. Keep container tightly closed.

P405 Store locked up.

P501 Dispose of contents/container to an approved waste disposal plant.

#### HMIS Classification

Health hazard: 2

Flammability: 1

Physical hazards: 0

#### NFPA Rating

Health Hazard: 2

Fire: 1

# SAFETY DATA SHEET

Cat# K163-100, CaspSELECT™ Caspase-3 Immunoassay Kit

SDS DATE: 13APR2015

Reactivity Hazard: 0

## Potential Health Effects

**Inhalation:** May be harmful if inhaled. Causes respiratory tract irritation.

**Skin:** Harmful if absorbed through skin. Causes skin irritation.

**Eyes:** Causes eye irritation.

**Ingestion:** Harmful if swallowed.

## DTT:

### Emergency Overview

**OSHA Hazards:** Target organ effect, Harmful by ingestion, Irritant

**Target Organs:** Central nervous system

**GHS Classification:** Acute toxicity, Oral (Category 4)

Skin irritation (Category 2)

Eye irritation (Category 2A)

**GHS Label elements, including precautionary statements**

**Pictogram:**



**Signal word:** Warning

**Hazard statement(s):** H302 Harmful if swallowed

H315 Causes skin irritation

H319 Causes serious eye irritation

H335: May cause respiratory irritation

**Precautionary statement(s):** P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

## HMIS Classification

**Health hazard:** 2

**Chronic Health Hazard:** \*

**Flammability:** 0

**Physical hazards:** 0

## NFPA Rating

**Health hazard:** 2

**Fire:** 0

**Reactivity Hazard:** 0

## Potential Health Effects

**Inhalation:** May be harmful if inhaled. Causes respiratory tract irritation.

**Skin:** May be harmful if absorbed through skin. Causes skin irritation.

**Eyes:** Causes eye irritation.

**Ingestion:** Harmful if swallowed.

## Bovine serum albumin (BSA):

**Emergency Overview:** Originates from USA where BSE has not been reported

**OSHA Hazards:** No known OSHA hazards

**GHS Classification:** Not a dangerous substance according to GHS

**GHS Label elements, including precautionary statements**

**Pictogram:** none

**Signal word:** none

**Hazard statement(s):** none

**Precautionary statement(s):** none

## HMIS Classification

**Health hazard:** 0

**Flammability:** 0

**Physical hazards:** 0

## NFPA Rating

**Health Hazard:** 0

**Fire:** 0

**Reactivity Hazard:** 0

## Potential Health Effects

**Inhalation:** May be harmful if inhaled. May cause respiratory tract irritation.

**Skin:** May be harmful if absorbed through skin. May cause skin irritation.

**Eyes:** May cause eye irritation.

**Ingestion:** May be harmful if swallowed.

## DMSO:

### Emergency Overview

**OSHA Hazards:** Combustible Liquid, Target Organ Effect

**Target Organs:** Eyes, Skin

**GHS Classification:** Flammable liquids (Category 4)

**GHS Label elements, including precautionary statements**

**Pictogram:** none

**Signal word:** Warning

**Hazard statement(s):** H227 Combustible liquid

# SAFETY DATA SHEET

Cat# K163-100, CaspSELECT™ Caspase-3 Immunoassay Kit

SDS DATE: 13APR2015

Precautionary statement(s): none

## HMIS Classification

Health hazard: 0

Chronic Health Hazard: \*

Flammability: 2

Physical hazards: 0

## NFPA Rating

Health hazard: 0

Fire: 2

Reactivity Hazard: 0

## Potential Health Effects

**Inhalation:** May be harmful if inhaled. May cause respiratory tract irritation.

**Skin:** May be harmful if absorbed through skin. May cause skin irritation.

**Eyes:** May cause eye irritation.

**Ingestion:** May be harmful if swallowed.

**Aggravated Medical Condition:** Avoid contact w/DMSO solutions containing toxic materials or materials with unknown toxicological properties. DMSO is readily absorbed through skin and may carry such materials into the body.

## SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS

Component	CAS Number	EC-No.	Molecular Weight	Chemical Formula
Triton X-100	9002-93-1	--	--	--
DTT	3483-12-3	222-468-7	154.25	C <sub>4</sub> H <sub>10</sub> O <sub>2</sub> S <sub>2</sub>
Glycerol	56-81-5	200-289-5	92.09	HOCH <sub>2</sub> CH(OH)CH <sub>2</sub> OH
Bovine Serum Albumin (BSA)	9048-46-8	232-936-2	--	--
DMSO	67-68-5	200-664-3	78.13	C <sub>2</sub> H <sub>6</sub> OS

## SECTION 4: FIRST AID MEASURES

**General advice:** Consult a physician. Show this safety data sheet to the doctor in attendance. Move out of dangerous area.

**If inhaled:** If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician.

**In case of skin contact:** Wash off with soap and plenty of water. Consult a physician.

**In case of eye contact:** Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician.

**If swallowed:** DO NOT induce vomiting. Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.

## SECTION 5: FIRE-FIGHTING MEASURES

### DMSO:

**Suitable extinguishing media:** For small (incipient) fires, use media such as "alcohol" foam, dry chemical, or carbon dioxide. For large fires, apply water from as far as possible. Use very large quantities (flooding) of water applied as a mist or spray; solid streams of water may be ineffective. Cool all affected containers with flooding quantities of water.

**Special protective equipment for firefighters:** Wear self-contained breathing apparatus for firefighting if necessary.

**Hazardous combustion products:** Hazardous combustion products formed under fire conditions – no data available.

**Further information:** Use water spray to cool unopened containers.

## SECTION 6: ACCIDENTAL RELEASE MEASURES

**Personal precautions:** Use personal protective equipment. Avoid dust formation. Avoid breathing vapors, mist, or gas. Ensure adequate ventilation. Evacuate personnel to safe areas. Avoid breathing dust.

**Environmental precautions:** Prevent further leakage or spillage if safe to do so. Do not let product enter drains. Discharge into the environment must be avoided.

**Methods for cleaning up:** Pick up and arrange disposal without creating dust. Sweep up and shovel. Keep in suitable, closed containers for disposal.

## SECTION 7: HANDLING AND STORAGE

### Precautions for safe handling

Avoid contact with skin and eyes. Avoid inhalation of vapor or mist. Use explosion-proof equipment. Keep away from sources of ignition (no smoking). Take measures to prevent the buildup of electrostatic charge.

### Conditions for safe storage

Keep container tightly closed in a dry and well-ventilated place.

Containers which are opened must be carefully resealed and kept upright to prevent leakage.

Recommended storage temperature: -20 °C

## SECTION 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

### Glycerol:

Components	CAS-No.	Value	Control parameters	Basis
------------	---------	-------	--------------------	-------

# SAFETY DATA SHEET

Cat# K163-100, CaspSELECT™ Caspase-3 Immunoassay Kit

SDS DATE: 13APR2015

Glycerol	56-81-5	TWA	10 mg/m <sup>3</sup>	USA. OSHA – TABLE Z-1 Limits for Air Contaminants – 1910.1000
		TWA	5 mg/m <sup>3</sup>	USA. OSHA – TABLE Z-1 Limits for Air Contaminants – 1910.1000
		TWA	10 mg/m <sup>3</sup>	USA. ACGIH Threshold Limit Values (TLV)
Remarks:	Upper Respiratory Tract Irritation			
		TWA	15 mg/m <sup>3</sup>	USA. Occupational Exposure Limits (OSHA) – Table Z-1 Limits for Air Contaminants
		TWA	5 mg/m <sup>3</sup>	USA. Occupational Exposure Limits (OSHA) – Table Z-1 Limits for Air Contaminants
	See Appendix D – Substances with no established RELs			

## DMSO:

Components	CAS-No.	Value	Control parameters	Basis
Dimethyl sulfoxide	67-68-5	TWA	250 ppm	USA. Workplace Environmental Exposure Levels (WEEL)

## Triton X-100, DTT, BSA:

Contains no substances with occupational exposure limit values.

### Personal protective equipment

#### Respiratory protection

Where risk assessment shows air-purifying respirators are appropriate use a full-face respirator with multi-purpose combination (US) or type ABEK (EN 14387) respirator cartridges as a backup to engineering controls. If the respirator is the sole means of protection, use a full-face supplied air respirator. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

#### Hand protection

Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.

#### Eye protection

Face shield and safety glasses. Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

#### Skin and body protection

Complete suit protecting against chemicals. The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

#### Hygiene measures

Handle in accordance with good industrial hygiene and safety practice. Wash hands before breaks and at the end of workday.

## SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

Property	Triton X-100	BSA	DTT	DMSO	Glycerol
Appearance:	Liquid	Lyophilized powder	Clear liquid	Clear liquid	Clear liquid
pH:	9.7	4.8-7.5	3.0-5.0	No data available	5.5-8
Water Solubility:	Soluble	Slightly soluble	Highly soluble	Completely miscible	Soluble
Other Solubility:	No data available	No data available	No data available	No data available	No data available
Boiling Point (°C)	>200 °C (>392 °F)	No data available	No data available	189 °C (372 °F)	182 °C (360 °F)
Melting Point (°C)	6 °C (43 °F)	No data available	No data available	16-19 °C (61-66 °F)	20 °C (68 °F)
Flash Point (°C)	251 °C (484 °F)	No data available	No data available	87 °C (189 °F)	160 °C (320 °F)
Ignition Temp. (°C)	No data available	No data available	No data available	301 °C (574 °F)	370 °C (698 °F)
Density:	1.070 g/ml	No data available	1.035 g/cm <sup>3</sup>	1.1 g/ml	1.25 g/ml

## SECTION 10: STABILITY AND REACTIVITY

Property	Triton X-100	BSA	DMSO	DTT	Glycerol
Chemical Stability:	Stable under recommended storage conditions				
Conditions to Avoid:	No data available	No data available	Heat, flames, sparks	No data available	No data available
Materials to Avoid:	Strong oxidizing agents, strong acids, strong bases	Strong oxidizing agents	Acid chlorides, phosphorus halides, strong acids, strong oxidizing agents, strong reducing agents	Bases, oxidizing agents, reducing agents, alkali metals	Strong bases, strong oxidizing agents
Hazardous decomposition products:	Carbon oxides	No data available	Carbon oxides, sulfur oxides	Carbon oxides, sulfur oxides, hydrogen sulfide gas	Carbon oxides

# SAFETY DATA SHEET

Cat# K163-100, CaspSELECT™ Caspase-3 Immunoassay Kit

SDS DATE: 13APR2015

## SECTION 11: TOXICOLOGICAL INFORMATION

### Triton X-100:

**Acute toxicity:** no data available

**Skin corrosion/irritation:** no data available

**Serious eye damage/eye irritation:** no data available

**Respiratory or skin sensitization:** no data available

**Germ cell mutagenicity:** no data available

### **Carcinogenicity:**

IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

ACGIH: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

NTP: No component of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

OSHA: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.

**Reproductive toxicity:** no data available

**Teratogenicity:** no data available

**Specific target organ toxicity – single exposure (GHS):** no data available

**Specific target organ toxicity – repeated exposure (GHS):** no data available

**Aspiration hazard:** no data available

### **Potential Health Effects**

**Inhalation:** May be harmful if inhaled. Causes respiratory tract irritation.

**Skin:** Harmful if absorbed through skin. Causes skin irritation.

**Eyes:** Causes eye irritation.

**Ingestion:** Harmful if swallowed.

**Signs and Symptoms of Exposure:** To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

**Synergistic effects:** no data available

**Additional information:** RTECS: not available

### DTT:

**Acute toxicity:** no data available

**Skin corrosion/irritation:** no data available

**Serious eye damage/eye irritation:** no data available

**Respiratory or skin sensitization:** no data available

**Germ cell mutagenicity:** no data available

### **Carcinogenicity:**

IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

ACGIH: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

NTP: No component of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

OSHA: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.

**Reproductive toxicity:** no data available

**Teratogenicity:** no data available

### **Potential Health Effects**

**Inhalation:** May be harmful if inhaled. Causes respiratory tract irritation.

**Skin:** May be harmful if absorbed through skin. Causes skin irritation.

**Eyes:** Causes eye irritation.

**Ingestion:** Harmful if swallowed.

**Signs and Symptoms of Exposure:** Exposure may cause central nervous system depression. To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

**Synergistic effects:** no data available

**Additional information:** RTECS: not available

### Glycerol:

**Acute toxicity:** LD50 Dermal – rabbit – >10,000 mg/kg

**Irritation and corrosion:** no data available

**Sensitization:** no data available

**Germ cell mutagenicity:** no data available

### **Carcinogenicity:**

IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

ACGIH: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

NTP: No component of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

OSHA: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.

**Reproductive toxicity:** no data available

# SAFETY DATA SHEET

Cat# K163-100, CaspSELECT™ Caspase-3 Immunoassay Kit

SDS DATE: 13APR2015

## Potential Health Effects

**Inhalation:** May be harmful if inhaled. May cause respiratory tract irritation.

**Skin:** May be harmful if absorbed through skin. May cause skin irritation.

**Eyes:** May cause eye irritation.

**Ingestion:** May be harmful if swallowed.

**Synergistic effects:** no data available

**Additional information:** RTECS: MA8050000

## DMSO:

**Acute toxicity:** LD50 Oral - rat - 14,500 mg/kg

LC50 Inhalation - rat - 4 h - 40250 ppm

LD50 Dermal - rabbit - > 5,000 mg/kg

**Skin corrosion/irritation:** no data available

**Serious eye damage/eye irritation:** no data available

**Respiratory/skin sensitization:** no data available

**Germ cell mutagenicity:** Genotoxicity in vitro - mouse – lymphocyte→ Cytogenetic analysis

Genotoxicity in vitro - mouse – lymphocyte→ Mutation in mammalian somatic cells.

Genotoxicity in vivo - rat – Intraperitoneal→ Cytogenetic analysis

Genotoxicity in vivo - mouse – Intraperitoneal→ DNA damage

**Carcinogenicity:** Carcinogenicity – rat – Oral→ Tumorigenic: Equivocal tumorigenic agent by RTECS criteria. Skin & Appendages: Other: Tumors.

Carcinogenicity – mouse – Oral→ Tumorigenic: Equivocal tumorigenic agent by RTECS criteria. Leukaemia Skin & Appendages: Other: Tumors.

IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

ACGIH: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

NTP: No component of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

OSHA: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.

**Reproductive toxicity:** Reproductive toxicity – rat – Intraperitoneal→ Effects on Fertility: Abortion.

Reproductive toxicity – rat – Intraperitoneal→ Effects on Fertility: Post-implantation mortality (e.g., dead and/or resorbed implants per total number of implants).

Reproductive toxicity – rat – Subcutaneous→ Effects on Fertility: Post-implantation mortality (e.g., dead and/or resorbed implants per total number of implants). Effects on Fertility: Litter size (e.g., # fetuses per litter; measured before birth).

Reproductive toxicity – mouse – Oral→ Effects on Fertility: Pre-implantation mortality (e.g., reduction in number of implants per female; total number of implants per corpora lutea). Effects on Embryo or Fetus: Fetotoxicity (except death, e.g., stunted fetus). Specific Developmental Abnormalities: Musculoskeletal system.

**Teratogenicity:** Developmental Toxicity – mouse – Intraperitoneal→ Effects on Embryo or Fetus: Fetotoxicity (except death, e.g., stunted fetus). Specific Developmental Abnormalities: Musculoskeletal system.

**Signs and Symptoms of Exposure:** Exposure via ingestion may cause nausea, fatigue, headache.

**Additional Information:** RTECS: PV6210000

---

## SECTION 12: ECOLOGICAL INFORMATION

---

### Triton X-100:

**Persistence and degradability:** no data available

**Toxicity:** no data available

**Bioaccumulative potential:** no data available

**Mobility in soil:** no data available

**PBT and vPvB assessment:** no data available

**Other adverse effects:** An environmental hazard cannot be excluded in the event of unprofessional handling or disposal.

Harmful to aquatic life with long lasting effects.

### DMSO:

**Persistence and degradability:** no data available

**Toxicity:** Toxicity to fish: LC50 – Pimephales promelas (fathead minnow) – 34,000 mg/l – 96 h

LC50 – Oncorhynchus mykiss (rainbow trout) – 35,000 mg/l – 96 h

Toxicity to daphnia and other aquatic invertebrates: EC50 – Daphnia pulex (Water flea) – 27,500 mg/l

Toxicity to algae: EC50 – Lepomis macrochirus (Bluegill) – >400,000 mg/l – 96 h

**Bioaccumulative potential:** no data available

**Mobility in soil:** no data available

**PBT and vPvB assessment:** no data available

**Other adverse effects:** no data available

---

## SECTION 13: DISPOSAL CONSIDERATIONS

---

**Product:** Observe all federal, state, and local environmental regulations. Contact a licensed professional waste disposal service to dispose of this material.

**Contaminated packaging:** Dispose of as unused product.

---

## SECTION 14: TRANSPORT INFORMATION

---

# SAFETY DATA SHEET

Cat# K163-100, CaspSELECT™ Caspase-3 Immunoassay Kit

SDS DATE: 13APR2015

Triton X-100, DTT, Glycerol, BSA:

DOT (US): Not dangerous goods.

IMDG: Not dangerous goods.

IATA: Not dangerous goods.

DMSO:

DOT (US): UN-Number: 1993, Class: CBL, Packing group: III; Proper shipping name: Combustible liquid, n.o.s. (Dimethyl sulfoxide); Marine pollutant: No; Poison Inhalation Hazard: No

IMDG: Not dangerous goods.

IATA: Not dangerous goods.

## SECTION 15: REGULATORY INFORMATION

**OSHA Hazards:** Triton X-100: Irritant, Harmful by ingestion

DTT: Target organ effect, Harmful by ingestion, Irritant

Glycerol: Target organ effect

DMSO: Combustible liquid, Target organ effect

**SARA 302 Components:** SARA 302: No chemical in this material are subject to the reporting requirements of SARA Title III, Section 302.

**SARA 313 Components:** SARA 313: This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title II, Section 313.

**SARA 311/312 Hazards:** Triton X-100: Acute Health Hazard

Glycerol: Chronic Health Hazard

DTT: Acute Health Hazard, Chronic Health Hazard

DMSO: Fire Hazard, Chronic Health Hazard

**Massachusetts Right To Know Components:** Glycerol, CAS-No. 56-81-5; Revision Date: 2007-03-01

**Pennsylvania Right To Know Components:** Triton X-100, CAS-No. 9002-93-1

DTT, CAS-No. 3483-12-3

Glycerol, CAS-No. 56-81-5; Revision Date: 2007-03-01

Serum albumin, CAS-No. 9048-46-8

Dimethyl sulfoxide CAS-No. 67-68-5; Revision Date: 2007-03-01

**New Jersey Right To Know Components:** Triton X-100, CAS-No. 9002-93-1

DTT, CAS-No. 3483-12-3

Glycerol, CAS-No. 56-81-5; Revision Date: 2007-03-01

Serum albumin, CAS-No. 9048-46-8

Dimethyl sulfoxide CAS-No. 67-68-5; Revision Date: 2007-03-01

**California Prop. 65 Components:** This product does not contain any chemicals known to State of California to cause cancer, birth defects, or any other reproductive harm.

## EU regulations

Component	Risk Phrases	Safety Phrases
Triton X-100	R22, R52/53	S22, S36/37/39, S45
DTT	R22, R36/37/38	S26, S36
Glycerol	--	S24/25
BSA	--	--
DMSO	R10, R36/37/38	S24/25, S36/37/39, S45

## SECTION 16: OTHER INFORMATION:

### OTHER INFORMATION:

### PREPARATION INFORMATION:

### DISCLAIMER:

The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. BioVision, Inc., shall not be held liable for any damage resulting from handling or from contact with the above product. See reverse side of invoice or packing slip for additional terms and conditions of sale.