BioVision

Recombinant Bovine FGF-2/FGF-basic

CATALOG #: 4040-10 10 μg

4040-50 50 μg 4040-1000 1 mg

LOT #:

SYNONYMS: HBGH-2, HBGF-2, Prostatropin, FGF-2, FGB-b

SOURCE: E. coli

PURITY: > 97 % by SDS-PAGE analyses and RP-HPLC

Endotoxin level is < 0.1 ng per μg.

MOL. WEIGHT: 17.25 kDa

FORM: Sterile filtered white lyophilized powder. Prepared fom a sterile solution

containing 1 % HSA.

STORAGE CONDITIONS: The lyophilized protein is best-stored desiccated at -20 °C.

RECONSTITUTION:

Centrifuge the vial prior to opening. Reconstitute in sterile water to a concentration not less than 0.1 mg/ml. Aliquot and store at -20°C. Avoid freeze-thaw cycles.

DESCRIPTION:

FGF-basic is a member of the fibroblast growth factor (FGF) family. FGF family members bind heparin and possess broad mitogenic and angiogenic activities. This protein has been implicated in diverse biological processes, such as limb and nervous system development, wound healing, and tumor growth. The heparin-binding growth factors are angiogenic agents in vivo and are potent mitogens for a variety of cell types in vitro. FGF-2 Bovine Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 155 amino acids and having a molecular mass of 17250 Dalton. The Fibroblast Growth Factor 2 is purified by proprietary chromatographic techniques. The sequence of the first five N-terminal amino acids was determined and was found to be Met-Ala-Ala-Gly-Ser.

BIOLOGICAL ACTIVITY:

The ED50, measured in a mitogenic assay using quiescent NR6R-3T3 fibroblasts was found to be <0.1 ng/ml, corresponding to a specific activity of 10 MUnits/mg.

FOR RESEARCH USE ONLY! Not to be used in human subjects!

RELATED PRODUCTS:

Growth Factors and Cytokines

Signal Transduction

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Monoclonal and Polyclonal Antibodies